

兔抗 SEPT1 多克隆抗体

中文名称：兔抗 SEPT1 多克隆抗体

英文名称： Anti-SEPT1 rabbit polyclonal antibody

别名： LARP; SEP1; DIFF6; PNUTL3

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： SEPT1

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis and the maintenance of cellular morphology. This gene encodes a protein that can form homo- and heterooligomeric filaments, and m

	<p>ay contribute to the formation of neurofibrillary tangles in Alzheimer's disease. Alternatively spliced transcript variants have been found but the full-length nature of these variants has not been determined.</p>
Applications:	WB
Name of antibody:	SEPT1
Immunogen:	Synthesized peptide derived from internal of human SEPT1.
Full name:	septin 1
Synonyms:	LARP; SEP1; DIFF6; PNUTL3
SwissProt:	Q8WYJ6
WB Predicted band size:	42 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

