

兔抗 PDS5A 多克隆抗体

中文名称：兔抗 PDS5A 多克隆抗体

英文名称： Anti-PDS5A rabbit polyclonal antibody

别名： PDS5 cohesin associated factor A; PIG54; SCC112; SCC-112

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PDS5A

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene binds to the cohesin complex and associates with chromatin through most of the cell cycle. The encoded protein may play a role in regulating sister chromatid cohesion during mitosis. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	PDS5A
Immunogen:	Fusion protein of human PDS5A
Full name:	PDS5 cohesin associated factor A
Synonyms:	PIG54; SCC112; SCC-112

SwissProt:	Q29RF7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer and Human lung cancer
IHC Recommend dilution:	50-200



