

兔抗 JAZF1 多克隆抗体

中文名称: 兔抗 JAZF1 多克隆抗体

英文名称: Anti-JAZF1 rabbit polyclonal antibody

别名: JAZF zinc finger 1; TIP27; ZNF802

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: JAZF1

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a nuclear protein with three C2H2-type zinc fingers, and functions as a transcriptional repressor. Chromosomal aberrations involving this gene are associated with endometrial stromal tumors. Alternatively spliced variants which encode different protein isoforms have been described; however, not all variants have been fully characterized.
Applications:	ELISA, IHC
Name of antibody:	JAZF1
Immunogen:	Synthetic peptide of human JAZF1
Full name:	JAZF zinc finger 1

Synonyms:	TIP27; ZNF802
SwissProt:	Q86VZ6
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	30-150

