

兔抗 JDP2 多克隆抗体

中文名称: 兔抗 JDP2 多克隆抗体

英文名称: Anti-JDP2 rabbit polyclonal antibody

别名: Jun dimerization protein 2; JUNDM2

储存: 冷冻 (-20℃)

抗原: JDP2

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Component of the AP-1 transcription factor that represses transcription mediated by the Jun family of proteins. Involved in a variety of transcriptional responses associated with AP-1 such as UV-induced apoptosis, cell differentiation, tumorigenesis and anti tumogeneris. Can also function as a repressor by recruiting histone deacetylase 3/HDAC3 to the promoter region of JUN. May control transcription via direct regulation of the modification of histones and the assembly of chromatin.
Applications:	ELISA, IHC
Name of antibody:	JDP2

Immunogen:	Synthetic peptide of human JDP2
Full name:	Jun dimerization protein 2
Synonyms:	JUNDM2
SwissProt:	Q8WYK2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human esophagus cancer
IHC Recommend dilution:	40-200



