

兔抗 SEPT3 多克隆抗体

中文名称: 兔抗 SEPT3 多克隆抗体

英文名称: Anti-SEPT3 rabbit polyclonal antibody

别名: SEP3; bK250D10.3

储存: 冷冻 (-20°C) 避光

抗原: SEPT3

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene belongs to the septin family of GTPases. Members of this family are required for cytokinesis. Expression is upregulated by retinoic acid in a human teratocarcinoma cell line. The specific function of this gene has not been determined. Alternative splicing of this gene results in two transcript variants encoding different isoforms.
Applications:	WB
Name of antibody:	SEPT3
Immunogen:	Synthesized peptide derived from internal of human SEPT3.

Full name:	septin 3
Synonyms:	SEP3; bK250D10.3
SwissProt:	Q9UH03
WB Predicted band size:	41 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution:	500-3000

