

兔抗 SEPT6 多克隆抗体

中文名称：兔抗 SEPT6 多克隆抗体

英文名称： Anti-SEPT6 rabbit polyclonal antibody

别名： SEP2; SEPT2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SEPT6

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis. One version of pediatric acute myeloid leukemia is the result of a reciprocal translocation between

	<p>chromosomes 11 and X, with the breakpoint associated with the genes encoding the mixed-lineage leukemia and septin 2 proteins. This gene encodes four transcript variants encoding three distinct isoforms. An additional transcript variant has been identified, but its biological validity has not been determined.</p>
Applications:	WB, IHC
Name of antibody:	SEPT6
Immunogen:	Synthesized peptide derived from internal of human SEPT6.
Full name:	septin 6
Synonyms:	SEP2; SEPT2
SwissProt:	Q14141
IHC positive control:	Human colon carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	50 kDa
WB Positive control:	Jurkat cells lysate

WB Recommended dilution: 500-3000



