

## 兔抗 JOSD1 多克隆抗体

中文名称：兔抗 JOSD1 多克隆抗体

英文名称：Anti-JOSD1 rabbit polyclonal antibody

别名：Josephin domain containing 1; dJ508I15.2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：JOSD1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Deubiquitinates monoubiquitinated probes (in vitro). When ubiquitinated, cleaves 'Lys-63'-linked and 'Lys-48'-linked poly-ubiquitin chains (in vitro), hence may act as a deubiquitinating enzyme. May increase macropinocytosis and suppress clathrin- and caveolae-mediated endocytosis. May enhance membrane dynamics and cell motility independently of its catalytic activity.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	JOSD1
<b>Immunogen:</b>	Fusion protein of human JOSD1
<b>Full name:</b>	Josephin domain containing 1
<b>Synonyms:</b>	dJ508I15.2

<b>SwissProt:</b>	Q15040
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	50-200

