

兔抗 SEPT14 多克隆抗体

中文名称：兔抗 SEPT14 多克隆抗体

英文名称： Anti-SEPT14 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

抗 原： SEPT14

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

SEPT14 is a member of the highly conserved septin family of GTP-binding cytoskeletal proteins implicated in membrane transport, apoptosis, cell polarity, cell cycle regulation, cytokinesis, and other cellular functions (Peterson et al., 2

	007 [PubMed 17922164]).
Applications:	WB
Name of antibody:	SEPT14
Immunogen:	Synthesized peptide derived from C-terminal of human SEPT14.
Full name:	septin 14
SwissProt:	Q6ZU15
WB Predicted band size:	50 kDa
WB Positive control:	NIH/3T3 cells lysate
WB Recommended dilution:	500-3000

