

兔抗 SEPTIN4 多克隆抗体

- 中文名称：兔抗 SEPTIN4 多克隆抗体
- 英文名称：Anti-SEPTIN4 rabbit polyclonal antibody
- 别名：septin 4; H5; ARTS; MART; SEP4; CE5B3; SEPT4; PNUT
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：SEPTIN4
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:

This gene is a member of the septin family of nucleotide binding proteins, originally described in yeast as cell division cycle regulatory proteins. Septins are highly conserved in yeast, Drosophila, and mouse, and appear to regulate cytoskeletal organization. Disruption of septin function disturbs cytokinesis and results in large multinucleate or polyploid cells. This gene is highly expressed in brain and heart. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. One of the isoforms (known as ARTS) is distinct; it is localized to the mitochondria, and

	has a role in apoptosis and cancer. [provided by RefSeq, Nov 2010]
Applications:	ELISA, IHC
Name of antibody:	SEPTIN4
Immunogen:	Synthetic peptide of human SEPTIN4
Full name:	septin 4
Synonyms:	H5; ARTS; MART; SEP4; CE5B3; SEPT4; PNUTL2; hucep-7; BRADEION; C17orf47; hCDCREL-2
SwissProt:	O43236
IHC positive control:	Human lung cancer and Human thyroid cancer
IHC Recommend dilution:	50-200

