

兔抗 SEPTIN9 多克隆抗体

中文名称：兔抗 SEPTIN9 多克隆抗体

英文名称： Anti-SEPTIN9 rabbit polyclonal antibody

别名： septin 9; MSF; MSF1; NAPB; SEPT9; SINT1; PNUTL4; S

抗原： SEPTIN9

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

IHC Recommend dilution:	25-100
IHC positive control:	Human cervical cancer and Human colon cancer
SwissProt:	Q9UHD8
Synonyms:	MSF; MSF1; NAPB; SEPT9; SINT1; PNUTL4; SeptD1; AF17q25
Full name:	septin 9
Immunogen:	Fusion protein of human SEPTIN9
Name of antibody:	SEPTIN9
Applications:	ELISA, IHC
Background:	This gene is a member of the septin family involved in cytokinesis and cell cycle control. This gene is a candidate

for the ovarian tumor suppressor gene. Mutations in this gene cause hereditary neuralgic amyotrophy, also known as neuritis with brachial predilection. A chromosomal translocation involving this gene on chromosome 17 and the MLL gene on chromosome 11 results in acute myelomonocytic leukemia. Multiple alternatively spliced transcript variants encoding different isoforms have been described.[provided by RefSeq, Mar 2009]

