

### 兔抗 SEPTIN10 多克隆抗体

中文名称：兔抗 SEPTIN10 多克隆抗体

英文名称： Anti-SEPTIN10 rabbit polyclonal antibody

别名： septin 10; SEPT10

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SEPTIN10

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the septin family of cytoskeletal proteins with GTPase activity. This protein localizes to the cytoplasm and nucleus and displays GTP-binding and GTPase activity. A pseudogene for this gene is located on chromosome 8. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2012]
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SEPTIN10
<b>Immunogen:</b>	Fusion protein of human SEPTIN10
<b>Full name:</b>	septin 10
<b>Synonyms:</b>	SEPT10
<b>SwissProt:</b>	Q9P0V9

<b>IHC positive control:</b>	Human cervical cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-300
<b>WB Predicted band size:</b>	53 kDa
<b>WB Positive control:</b>	Human fetal liver tissue lysate
<b>WB Recommended dilution:</b>	200-1000



