

兔抗 WDHD1 多克隆抗体

中文名称：兔抗 WDHD1 多克隆抗体

英文名称：Anti-WDHD1 rabbit polyclonal antibody

级别：AND1; CTF4; AND-1; CHTF4

储存：冷冻（-20℃）

抗原：WDHD1

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene contains multiple N-terminal WD40 domains and a C-terminal high mobility group (HMG) box. WD40 domains are found in a variety of eukaryotic proteins and may function as adaptor/regulatory modules in signal transduction, pre-mRNA processing and cytoskeleton assembly. HMG boxes are found in many eukaryotic proteins involved in chromatin assembly, transcription and replication. Alternative splicing results in two transcript variants encoding different isoforms.
Applications:	ELISA, IHC

Name of antibody:	WDHD1
Immunogen:	Synthetic peptide of human WDHD1
Full name:	WD repeat and HMG-box DNA binding protein 1
Synonyms:	AND1; CTF4; AND-1; CHTF4
SwissProt:	O75717
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	25-100



