

## 兔抗 WDR1 多克隆抗体

中文名称：兔抗 WDR1 多克隆抗体

英文名称： Anti-WDR1 rabbit polyclonal antibody

别名： WD repeat domain 1; AIP1; NORI-1; HEL-S-52

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： WDR1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein containing 9 WD repeats. WD repeats are approximately 30- to 40-amino acid domains containing several conserved residues, mostly including a trp-asp at the C-terminal end. WD domains are involved in protein-protein interactions. The encoded protein may help induce the disassembly of actin filaments. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	WDR1
<b>Immunogen:</b>	Fusion protein of human WDR1

<b>Full name:</b>	WD repeat domain 1
<b>Synonyms:</b>	AIP1; NORI-1; HEL-S-52
<b>SwissProt:</b>	O75083
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human tonsil
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	66 kDa
<b>WB Positive control:</b>	A172 cell lysate
<b>WB Recommended dilution:</b>	200-1000



