

兔抗 PEBP1 多克隆抗体

中文名称：兔抗 PEBP1 多克隆抗体

英文名称： Anti-PEBP1 rabbit polyclonal antibody

别名： PBP; HCNP; PEBP; RKIP; HCNPpp; PEBP-1

储存： 冷冻 (-20℃)

抗原： PEBP1

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Phosphatidylethanolamine-binding protein 1 is a protein that in humans is encoded by the PEBP1 gene. Binds ATP, opioids and phosphatidylethanolamine. Has lower affinity for phosphatidylinositol and phosphatidylcholine. Serine protease inhibitor which inhibits thrombin, neurosin and chymotrypsin but not trypsin, tissue type plasminogen activator and elastase By similarity. Inhibits the kinase activity of RAF1 by inhibiting its activation and by dissociating the RAF1/MEK complex and acting as a competitive inhibitor of MEK phosphorylation.
Applications:	ELISA, WB, IHC

Name of antibody:	PEBP1
Immunogen:	Fusion protein of human PEBP1
Full name:	phosphatidylethanolamine binding protein 1
Synonyms :	PBP; HCNP; PEBP; RKIP; HCNPPP; PEBP-1
SwissProt:	P30086
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer and human lung cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	21 kDa
WB Positive control:	A549, K562, hela and hepG2 cells
WB Recommended dilution:	500-2000



