

兔抗 SERP1 多克隆抗体

中文名称：兔抗 SERP1 多克隆抗体

英文名称： Anti-SERP1 rabbit polyclonal antibody

别名： RAMP4

抗原： SERP1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene interacts with target proteins during their translocation into the lumen of the endoplasmic reticulum. Protects unfolded target proteins against degradation during ER stress. May facilitate glycosylation of target proteins after termination of ER stress. May modulate the use of N-glycosylation sites on target proteins.
Applications:	ELISA, IHC
Name of antibody:	SERP1
Immunogen:	Fusion protein of human SERP1
Full name:	stress-associated endoplasmic reticulum protein 1

Synonyms :	RAMP4
SwissProt:	Q9Y6X1
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human colon cancer and human lung cancer
IHC Recommend dilution:	25-100

