

兔抗 SERPINA3 多克隆抗体

中文名称：兔抗 SERPINA3 多克隆抗体

英文名称：Anti-SERPINA3 rabbit polyclonal antibody

别名：ACT; AACT; GIG24; GIG25

相关类别：一抗

储存：冷冻（-20℃）

抗原：SERPINA3

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a plasma protease inhibitor and member of the serine protease inhibitor class. Polymorphisms in this protein appear to be tissue specific and influence protease targeting. Variations in this protein's sequence have been implicated in Alzheimer's disease, and deficiency of this protein has been associated with liver disease. Mutations have been identified in patients with Parkinson disease and chronic obstructive pulmonary disease.
Applications:	ELISA, IHC
Name of antibody:	SERPINA3

Immunogen:	Fusion protein of human SERPINA3
Full name:	serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 3
Synonyms :	ACT; AACT; GIG24; GIG25
SwissProt:	P01011
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	50-200

