

## 兔抗 WDR6 多克隆抗体

中文名称：兔抗 WDR6 多克隆抗体

英文名称：Anti-WDR6 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：WDR6

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-aspartic acid (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. The encoded protein interacts with serine/threonine kinase 11, and is implicated in cell growth arrest. Alternative splicing results in multiple transcript variants encoding different isoforms.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	WDR6
<b>Immunogen:</b>	Fusion protein of human WDR6

<b>Full name:</b>	WD repeat domain 6
<b>SwissProt:</b>	Q9NNW5
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer;
<b>IHC Recommend dilution:</b>	30-150

