

兔抗 SERPINA9 多克隆抗体

中文名称：兔抗 SERPINA9 多克隆抗体

英文名称：Anti-SERPINA9 rabbit polyclonal antibody

别名：serpin family A member 9; GCET1; SERPINA11; SERPINA11b

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SERPINA9

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Protease inhibitor that inhibits trypsin and trypsin-like serine proteases (in vitro). Inhibits plasmin and thrombin with lower efficiency (in vitro).
Applications:	ELISA, WB, IHC
Name of antibody:	SERPINA9
Immunogen:	Synthetic peptide of human SERPINA9
Full name:	serpin family A member 9
Synonyms:	GCET1; SERPINA11; SERPINA11b
SwissProt:	Q86WD7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer

IHC Recommend dilution:	20-100
WB Predicted band size:	47 kDa
WB Positive control:	Raji cell lysate
WB Recommended dilution:	200-1000



