

## 兔抗 PEF1 多克隆抗体

中文名称：兔抗 PEF1 多克隆抗体

英文名称：Anti-PEF1 rabbit polyclonal antibody

别名：ABP32; PEF1A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PEF1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a calcium-binding protein belonging to the penta-EF-hand protein family. The encoded protein has been shown to form a heterodimer with the programmed cell death 6 gene product and may modulate its function in Ca <sup>2+</sup> signaling. Alternative splicing results in multiple transcript variants and a pseudogene has been identified on chromosome 1.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PEF1
<b>Immunogen:</b>	Fusion protein of human PEF1
<b>Full name:</b>	penta-EF-hand domain containing 1

<b>Synonyms :</b>	ABP32; PEF1A
<b>SwissProt:</b>	Q9UBV8
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	30 kDa
<b>WB Positive control:</b>	Jurkat and hela cells, 293T cells and human thyroid cancer tissue
<b>WB Recommended dilution:</b>	500-2000



