

兔抗 SERPINA11 多克隆抗体

中文名称：兔抗 SERPINA11 多克隆抗体

英文名称： Anti-SERPINA11 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： SERPINA11

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein belongs to the serpin family. SERPINA11 has serine-type endopeptidase inhibitor activity, negative regulation of endopeptidase activity, regulation of proteolysis, serine-type endopeptidase inhibitor activity.
Applications:	ELISA, IHC
Name of antibody:	SERPINA11
Immunogen:	Synthetic peptide of human SERPINA11
Full name:	serpin peptidase inhibitor, clade A (alpha-1 antitrypsin), member 11
SwissProt:	Q86U17
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human breast cancer

IHC Recommend dilution:

25-100

