

兔抗 WDR18 多克隆抗体

中文名称: 兔抗 WDR18 多克隆抗体

英文名称: Anti-WDR18 rabbit polyclonal antibody

别名: Ipi3; R32184_1

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: WDR18

反应种属: Human, Mouse

标记物: Unconjugate

相关类别: 一抗

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. Members of this family are involved in a variety of cellular processes, including cell cycle progression, signal transduction, apoptosis, and gene regulation.
Applications:	ELISA, IHC
Name of antibody:	WDR18

Immunogen:	Fusion protein of human WDR18
Full name:	WD repeat domain 18
Synonyms :	Ipi3; R32184_1
SwissProt:	Q9BV38
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human breast cancer
IHC Recommend dilution:	25-100



