

兔抗 SERPINA12 多克隆抗体

中文名称：兔抗 SERPINA12 多克隆抗体

英文名称： Anti-SERPINA12 rabbit polyclonal antibody

别名： OL-64

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SERPINA12

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Adipokine that modulates insulin action by specifically inhibiting its target protease KLK7 in white adipose tissues.
Applications:	ELISA, IHC
Name of antibody:	SERPINA12
Immunogen:	Fusion protein of human SERPINA12
Full name:	serpin peptidase inhibitor, clade A (alpha-1 antiprotease, antitrypsin), member 12
Synonyms :	OL-64
SwissProt:	Q8IW75
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human lung cancer and human liver cancer
IHC Recommend dilution:	25-100

