

兔抗 JPT2 多克隆抗体

中文名称：兔抗 JPT2 多克隆抗体

英文名称：Anti-JPT2 rabbit polyclonal antibody

别名：Jupiter microtubule associated homolog 2; L11; HN1L; C16orf34

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：JPT2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	JPT2
Immunogen:	Fusion protein of human JPT2
Full name:	Jupiter microtubule associated homolog 2
Synonyms:	L11; HN1L; C16orf34
SwissProt:	Q9H910
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	30-150

