

兔抗 JRKL 多克隆抗体

中文名称：兔抗 JRKL 多克隆抗体

英文名称：Anti-JRKL rabbit polyclonal antibody

别名：JRK like; HHMJG

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：JRKL

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this gene has not yet been defined, however, the encoded protein shares similarity with the human (41% identical) and mouse (34% identical) jerky gene products. This protein may act as a nuclear regulatory protein. Alternatively spliced transcript variants have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	JRKL
Immunogen:	Synthetic peptide of human JRKL
Full name:	JRK like
Synonyms:	HHMJG
SwissProt:	Q9Y4A0
ELISA Recommended dilution:	5000-10000

WB Predicted band size:	60 kDa
WB Positive control:	Juakat and K562 cell lysates
WB Recommended dilution:	1000-5000

