

## 兔抗 SERPINB2 多克隆抗体

中文名称：兔抗 SERPINB2 多克隆抗体

英文名称： Anti-SERPINB2 rabbit polyclonal antibody

别名： PAI; PAI2; PAI-2; PLANH2; HsT1201

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SERPINB2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Plasminogen activator inhibitor-2 (placental PAI) is a coagulation factor that inactivates tPA and urokinase. It is present in most cells, especially monocytes/macrophages. PAI-2 exists in two forms, a 60-kDa extracellular glycosylated form and a 43-kDa intracellular form. It is present only at detectable quantities in blood during pregnancy, as it is produced by the placenta, and may explain partially the increased rate of thrombosis during pregnancy.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SERPINB2

<b>Immunogen:</b>	Fusion protein of human SERPINB2
<b>Full name:</b>	serpin peptidase inhibitor, clade B (ovalbumin), member 2
<b>Synonyms :</b>	PAI; PAI2; PAI-2; PLANH2; HsT1201
<b>SwissProt:</b>	P05120
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human ovarian cancer
<b>IHC Recommend dilution:</b>	25-100

