

兔抗 SERPINB4 多克隆抗体

中文名称：兔抗 SERPINB4 多克隆抗体

英文名称： Anti-SERPINB4 rabbit polyclonal antibody

别名： serpin family B member 4; PI11; SCCA1; SCCA2; LEUPIN; SCCA-2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SERPINB4

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the serpin family of serine protease inhibitors. The encoded protein is highly expressed in many tumor cells and can inactivate granzyme M, an enzyme that kills tumor cells. This protein, along with serpin B3, can be processed into smaller fragments that aggregate to form an autoantigen in psoriasis, probably by causing chronic inflammation.
Applications:	ELISA, IHC
Name of antibody:	SERPINB4
Immunogen:	Fusion protein of human SERPINB4

Full name:	serpin family B member 4
Synonyms:	PI11; SCCA1; SCCA2; LEUPIN; SCCA-2
SwissProt:	P48594
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-100

