

## 兔抗 JUN(Ab-243) 多克隆抗体

- 中文名称: 兔抗 JUN(Ab-243) 多克隆抗体
- 英文名称: Anti-JUN(Ab-243) rabbit polyclonal antibody
- 别名: AP1; AP-1; c-Jun
- 相关类别: 一抗
- 储存: 冷冻 (-20℃) 避光
- 宿主: Rabbit
- 抗原: JUN(Ab-243)
- 反应种属: Human Mouse Rat
- 标记物: Unconjugate
- 克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	JUN(Ab-243)
<b>Immunogen:</b>	Synthesized non-phosphopeptide derived from human c-Jun around the phosphorylation site of serine 243.
<b>Full name:</b>	jun proto-oncogene
<b>Synonyms :</b>	AP1; AP-1; c-Jun

<b>SwissProt:</b>	P05412
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	36 kDa
<b>WB Positive control:</b>	Hela cells lysate
<b>WB Recommended dilution:</b>	500-3000



