

兔抗 PEX10 多克隆抗体

中文名称: 兔抗 PEX10 多克隆抗体

英文名称: Anti-PEX10 rabbit polyclonal antibody

别名: NALD; PBD6A; PBD6B; RNF69

抗原: PEX10

储存: 冷冻 (-20°C) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene encodes a protein involved in import of peroxisomal matrix proteins. This protein localizes to the peroxisomal membrane. Mutations in this gene result in phenotypes within the Zellweger spectrum of peroxisomal biogenesis disorders, ranging from neonatal adrenoleukodystrophy to Zellweger syndrome. Alternative splicing results in two transcript variants encoding different isoforms.
Applications:	WB, IHC
Name of antibody:	PEX10
Immunogen:	Synthesized peptide derived from internal of human PEX10.
Full name:	peroxisomal biogenesis factor 10

Synonyms :	NALD; PBD6A; PBD6B; RNF69
SwissProt:	O60683
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	37 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000



