

## 兔抗 SERPINB6 多克隆抗体

中文名称：兔抗 SERPINB6 多克隆抗体

英文名称： Anti-SERPINB6 rabbit ployclonal antibody

别名： CAP; PI6; PTI; PI-6; SPI3; DFNB91; MSTP057

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SERPINB6

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

The protein encoded by this gene is a member of the serpin (serine proteinase inhibitor) superfamily, and ovalbumin(ov)-serpin subfamily. It was originally discovered as a placental thrombin inhibitor. The mouse homolog was found to be expressed in the hair cells of the inner ear. Mutations in this gene are associated with nonsyndromic progressive hearing loss, suggesting that this serpin plays an important role in the inner ear in the protection against leakage of lysosomal content during stress, and that loss of this protection results in cell

	death and sensorineural hearing loss. Alternatively spliced transcript variants have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SERPINB6
<b>Immunogen:</b>	Fusion protein of human SERPINB6
<b>Full name:</b>	serpin peptidase inhibitor, clade B (ovalbumin), member 6
<b>Synonyms :</b>	CAP; PI6; PTI; PI-6; SPI3; DFNB91; MSTP057
<b>SwissProt:</b>	P35237
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human breast cancer
<b>IHC Recommend dilution:</b>	25-100

