

## 兔抗 SERPINB7 多克隆抗体

中文名称：兔抗 SERPINB7 多克隆抗体

英文名称：Anti-SERPINB7 rabbit polyclonal antibody

别名：serpin family B member 7; PPKN; TP55; MEGSIN

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SERPINB7

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a family of proteins which function as protease inhibitors. Expression of this gene is upregulated in IgA nephropathy and mutations have been found to cause palmoplantar keratoderma, Nagashima type. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SERPINB7
<b>Immunogen:</b>	Synthetic peptide of human SERPINB7
<b>Full name:</b>	serpin family B member 7
<b>Synonyms:</b>	PPKN; TP55; MEGSIN

<b>SwissProt:</b>	O75635
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	43 kDa
<b>WB Positive control:</b>	HL60 cell lysate
<b>WB Recommended dilution:</b>	500-2000



