

兔抗 JUN (Phospho-Ser73)多克隆抗体

中文名称:	兔抗 JUN (Phospho-Ser73)多克隆抗体
英文名称:	Anti-JUN (Phospho-Ser73) rabbit polyclonal antibody
别名:	AH119; AP1; Jun A; c-Jun; p39
相关类别:	一抗
储存:	冷冻 (-20℃) 避光
宿主:	Rabbit
抗原:	JUN (Phospho-Ser73)
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.
Applications:	WB, IHC
Name of antibody:	JUN (Phospho-Ser73)
Immunogen:	Synthetic peptide of human JUN (Phospho-Ser73)
Full name:	jun proto-oncogene (Phospho-Ser73)
Synonyms:	AH119; AP1; Jun A; c-Jun; p39

SwissProt:	P05412
IHC positive control:	Human breast carcinoma
IHC Recommend dilution:	50-100
WB Predicted band size:	36 kDa
WB Positive control:	PC14 cells treated with EGF
WB Recommended dilution:	500-1000



