

兔抗 PEX2 多克隆抗体

中文名称: 兔抗 PEX2 多克隆抗体

英文名称: Anti-PEX2 rabbit polyclonal antibody

别名: PAF1; PMP3; ZWS3; PBD5A; PBD5B; PMP35; PXMP3; RNF72

储存: 冷冻 (-20°C) 避光

抗原: PEX2

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene encodes an integral peroxisomal membrane protein required for peroxisome biogenesis. The protein is thought to be involved in peroxisomal matrix protein import. Mutations in this gene result in one form of Zellweger syndrome and infantile Refsum disease. Alternative splicing results in multiple transcript variants encoding the same protein.
Applications:	WB, IHC
Name of antibody:	PEX2
Immunogen:	Synthesized peptide derived from N-terminal of human PXP3.
Full name:	peroxisomal biogenesis factor 2

Synonyms :	PAF1; PMP3; ZWS3; PBD5A; PBD5B; PMP35; PXMP3; RNF72
SwissProt:	P28328
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	35 kDa
WB Positive control:	A549 cells lysate
WB Recommended dilution:	500-3000



