

兔抗 SERPINB8 多克隆抗体

中文名称：兔抗 SERPINB8 多克隆抗体

英文名称： Anti-SERPINB8 rabbit ployclonal antibody

别名： PI8; CAP2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SERPINB8

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The superfamily of high molecular weight serine protease inhibitors (serpins) regulate a diverse set of intracellular and extracellular processes such as complement activation, fibrinolysis, coagulation, cellular differentiation, tumor suppression, apoptosis, and cell migration. Serpins are characterized by well-conserved a tertiary structure that consists of 3 beta sheets and 8 or 9 alpha helices.
Applications:	ELISA, WB, IHC
Name of antibody:	SERPINB8
Immunogen:	Full length fusion protein

Full name:	serpin peptidase inhibitor, clade B (ovalbumin), member 8
Synonyms :	PI8; CAP2
SwissProt:	P50452
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	43 kDa
WB Positive control:	Human normal liver tissue
WB Recommended dilution:	200-1000



