

## 兔抗 SERPINB9 多克隆抗体

中文名称：兔抗 SERPINB9 多克隆抗体

英文名称： Anti-SERPINB9 rabbit polyclonal antibody

别名： PI9; CAP3; PI-9; CAP-3

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SERPINB9

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

**Background:**

This gene encodes a member of the serine protease inhibitor family which are also known as serpins. The encoded protein belongs to a subfamily of intracellular serpins. This protein i

	<p>inhibits the activity of the effector molecule granzyme B. Overexpression of this protein may prevent cytotoxic T-lymphocytes from eliminating certain tumor cells. A pseudogene of this gene is found on chromosome 6.</p>
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	SERPINB9
<b>Immunogen:</b>	Synthesized peptide derived from internal of human SERPINB9.
<b>Full name:</b>	serpin peptidase inhibitor, clade B (ovalbumin), member 9
<b>Synonyms :</b>	PI9; CAP3; PI-9; CAP-3
<b>SwissProt:</b>	P50453
<b>IHC positive control:</b>	Human colon carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	42 kDa
<b>WB Positive control:</b>	Jurkat cells lysate
<b>WB Recommended dilution:</b>	500-3000

