

## 兔抗 JUN 多克隆抗体

中文名称： 兔抗 JUN 多克隆抗体

英文名称： Anti-JUN rabbit polyclonal antibody

别名： AP1; AP-1; c-Jun

相关类别： 一抗

抗原： JUN

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat


标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intro nless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	JUN
<b>Immunogen:</b>	Fusion protein of human JUN
<b>Full name:</b>	jun proto-oncogene

<b>Synonyms :</b>	AP1; AP-1; c-Jun
<b>SwissProt:</b>	P05412
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human lung cancer and Human colon cancer
<b>IHC Recommend dilution:</b>	25-100



**mlbio** 梅联生物  
Good elisakit producers

