

## 兔抗 PFDN1 多克隆抗体

中文名称：兔抗 PFDN1 多克隆抗体

英文名称： Anti-PFDN1 rabbit polyclonal antibody

别名： PDF; PFD1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： PFDN1

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a member of the prefoldin beta subunit family. The encoded protein is one of six subunits of prefoldin, a molecular chaperone complex that binds and stabilizes newly synthesized polypeptides, thereby allowing them to fold correctly. The complex, consisting of two alpha and four beta subunits, forms a double beta barrel assembly with six protruding coiled-coils.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	PFDN1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human PFDN1.

<b>Full name:</b>	prefoldin subunit 1
<b>Synonyms :</b>	PDF; PFD1
<b>SwissProt:</b>	O60925
<b>WB Predicted band size:</b>	13 kDa
<b>WB Positive control:</b>	RAW264.7 cells lysate
<b>WB Recommended dilution:</b>	500-3000

