

兔抗 PFDN2 多克隆抗体

中文名称：兔抗 PFDN2 多克隆抗体

英文名称： Anti-PFDN2 rabbit polyclonal antibody

别名： prefoldin subunit 2; PFD2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PFDN2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the prefoldin beta subunit family. The encoded protein is one of six subunits of prefoldin, a molecular chaperone complex that binds and stabilizes newly synthesized polypeptides, thereby allowing them to fold correctly. The complex, consisting of two alpha and four beta subunits, forms a double beta barrel assembly with six protruding coiled-coils.
Applications:	ELISA, IHC
Name of antibody:	PFDN2
Immunogen:	Fusion protein of human PFDN2
Full name:	prefoldin subunit 2

Synonyms:	PFD2
SwissProt:	Q9UHV9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human breast cancer
IHC Recommend dilution:	25-100



