

## 兔抗 PFDN6 多克隆抗体

中文名称：兔抗 PFDN6 多克隆抗体

英文名称：Anti-PFDN6 rabbit polyclonal antibody

别名：HKE2; KE-2; PFD6; H2-KE2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PFDN6

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	PFDN6 is a subunit of the heteromeric prefoldin complex that chaperones nascent actin (see MIM 102560) and alpha- and beta-tubulin (see MIM 602529 and MIM 191130, respectively) chains pending their transfer to the cytosolic chaperonin containing TCP1 (MIM 186980) (CCT) complex (Hansen et al., 1999 [PubMed 10209023]).
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PFDN6
<b>Immunogen:</b>	Synthetic peptide of human PFDN6
<b>Full name:</b>	prefoldin subunit 6
<b>Synonyms :</b>	HKE2; KE-2; PFD6; H2-KE2

<b>SwissProt:</b>	O15212
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human lung cancer and Human colon cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	14 kDa
<b>WB Positive control:</b>	Mouse heart tissue
<b>WB Recommended dilution:</b>	500-2000



