

兔抗 WDR47 多克隆抗体

中文名称：兔抗 WDR47 多克隆抗体

英文名称：Anti-WDR47 rabbit polyclonal antibody

别名：WD repeat domain 47

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

抗原：WDR47

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	WDR47
Immunogen:	Fusion protein of human WDR47
Full name:	WD repeat domain 47
SwissProt:	O94967
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	30-150

