

兔抗 SERPINF1 多克隆抗体

中文名称：兔抗 SERPINF1 多克隆抗体

英文名称： Anti-SERPINF1 rabbit polyclonal antibody

别名： OI6, OI12, PEDF, EPC-1, PIG35

抗原： SERPINF1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the serpin family, although it does not display the serine protease inhibitory activity shown by many of the other serpin family members. The encoded protein is secreted and strongly inhibits angiogenesis. In addition, this protein is a neurotrophic factor involved in neuronal differentiation in retinoblastoma cells.
Applications:	ELISA, WB, IHC
Name of antibody:	SERPINF1
Immunogen:	Fusion protein of human SERPINF1
Full name:	serpin peptidase inhibitor, clade F (alpha-2 antiplasmin,

	pigment epithelium derived factor), member 1
Synonyms :	OI6, OI12, PEDF, EPC-1, PIG35
SwissProt:	P36955
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human tonsil
IHC Recommend dilution:	15-50
WB Predicted band size:	46 kDa
WB Positive control:	K562 cells
WB Recommended dilution:	200-1000



