

## 兔抗 SERPINF2 多克隆抗体

中文名称：兔抗 SERPINF2 多克隆抗体

英文名称：Anti-SERPINF2 rabbit polyclonal antibody

别名：AAP; API; PLI; A2AP; ALPHA-2-PI

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SERPINF2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the serpin family of serine protease inhibitors. The protein is a major inhibitor of plasmin, which degrades fibrin and various other proteins. Consequently, the proper function of this gene has a major role in regulating the blood clotting pathway. Mutations in this gene result in alpha-2-plasmin inhibitor deficiency, which is characterized by severe hemorrhagic diathesis. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SERPINF2

<b>Immunogen:</b>	Synthetic peptide of human SERPINF2
<b>Full name:</b>	serpin peptidase inhibitor, clade F (alpha-2 antiplasmin, pigment epithelium derived factor), member 2
<b>Synonyms :</b>	AAP; API; PLI; A2AP; ALPHA-2-PI
<b>SwissProt:</b>	P08697
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	55 kDa
<b>WB Positive control:</b>	Human normal stomach tissue
<b>WB Recommended dilution:</b>	500-2000



