

兔抗 KANK3 多克隆抗体

- 中文名称：兔抗 KANK3 多克隆抗体
- 英文名称：Anti-KANK3 rabbit polyclonal antibody
- 别名：KN motif and ankyrin repeat domains 3; ANKRD47
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：KANK3
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in the control of cytoskeleton formation by regulating actin polymerization.
Applications:	ELISA, IHC
Name of antibody:	KANK3
Immunogen:	Synthetic peptide of human KANK3
Full name:	KN motif and ankyrin repeat domains 3
Synonyms:	ANKRD47
SwissProt:	Q6NY19
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer

IHC Recommend dilution:	50-300
-------------------------	--------

