

兔抗 SERPINI1 多克隆抗体

中文名称：兔抗 SERPINI1 多克隆抗体

英文名称：Anti-SERPINI1 rabbit polyclonal antibody

别名：PI12; SERPINI1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：SERPINI1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Serine protease inhibitor that inhibits plasminogen activators and plasmin but not thrombin. May be involved in the formation or reorganization of synaptic connections as well as for synaptic plasticity in the adult nervous system. May protect neurons from cell damage by tissue-type plasminogen activator.
Applications:	WB
Name of antibody:	SERPINI1
Immunogen:	Synthetic peptide of human SERPINI1
Full name:	serpin peptidase inhibitor, clade I (neuroserpin), member 1
Synonyms :	PI12; SERPINI1

SwissProt:	Q99574
WB Predicted band size:	46 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	500-1000

