

兔抗 SETD3 多克隆抗体

中文名称：兔抗 SETD3 多克隆抗体

英文名称：Anti-SETD3 rabbit polyclonal antibody

别名：C14orf154

相关类别：一抗

抗原：SETD3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Histone methyltransferase that methylates 'Lys-4' and 'Lys-36' of histone H3 (H3K4me and H3K36me). Acts as a transcriptional activator. Plays an important role in the transcriptional regulation of muscle cell differentiation via interaction with MYOD1.
Applications:	ELISA, WB, IHC
Name of antibody:	SETD3
Immunogen:	Fusion protein of human SETD3
Full name:	SET domain containing 3
Synonyms :	C14orf154
SwissProt:	Q86TU7

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human ovarian cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	67 kDa
WB Positive control:	Hela , 293T and HEPG2 cell
WB Recommended dilution:	200-1000



