

兔抗 SETDB1 多克隆抗体

中文名称：兔抗 SETDB1 多克隆抗体

英文名称： Anti-SETDB1 rabbit polyclonal antibody

别名： SET domain bifurcated 1; ESET; KG1T; KMT1E; TDRD21; H3-K9-HMTase4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SETDB1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a histone methyltransferase which regulates histone methylation, gene silencing, and transcriptional repression. This gene has been identified as a target for treatment in Huntington Disease, given that gene silencing and transcription dysfunction likely play a role in the disease pathogenesis. Alternatively spliced transcript variants of this gene have been described.
Applications:	ELISA, IHC
Name of antibody:	SETDB1
Immunogen:	Fusion protein of human SETDB1

Full name:	SET domain bifurcated 1
Synonyms:	ESET; KG1T; KMT1E; TDRD21; H3-K9-HMTase4
SwissProt:	Q15047
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer Human thyroid cancer
IHC Recommend dilution:	50-300



