

## 兔抗 WEE1 多克隆抗体

中文名称：兔抗 WEE1 多克隆抗体

英文名称：Anti-WEE1 rabbit polyclonal antibody

别名：WEE1 G2 checkpoint kinase WEE1A; WEE1hu

抗原：WEE1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a nuclear protein, which is a tyrosine kinase belonging to the Ser/Thr family of protein kinases. This protein catalyzes the inhibitory tyrosine phosphorylation of CDC2/cyclin B kinase, and appears to coordinate the transition between DNA replication and mitosis by protecting the nucleus from cytoplasmically activated CDC2 kinase.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	WEE1
<b>Immunogen:</b>	Synthetic peptide of human WEE1
<b>Full name:</b>	WEE1 G2 checkpoint kinase

<b>Synonyms:</b>	WEE1A; WEE1hu
<b>SwissProt:</b>	P30291
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer; Human ovarian cancer
<b>IHC Recommend dilution:</b>	25-100



