

兔抗 SF1(Ab-82) 多克隆抗体

中文名称：兔抗 SF1(Ab-82) 多克隆抗体

英文名称： Anti-SF1(Ab-82) rabbit polyclonal antibody

别名： BBP; MBBP; ZFM1; ZNF162; D11S636; ZCCHC25

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SF1(Ab-82)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a nuclear pre-mRNA splicing factor. The encoded protein specifically recognizes the intron branch point sequence and is required for the early stages of spliceosome assembly. Alternate splicing results in

	multiple transcript variants.
Applications:	WB, IHC, IF
Name of antibody:	SF1(Ab-82)
Immunogen:	Synthesized non-phosphopeptide derived from human SF1 around the phosphorylation site of serine 82 (S-P-S(p)-P-E).
Full name:	splicing factor 1
Synonyms :	BBP; MBBP; ZFM1; ZNF162; D11S636; ZCCHC25
SwissProt:	Q15637
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	68 kDa
WB Positive control:	COLO205 cells lysate
WB Recommended dilution:	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500



