

## 兔抗 PGBD2 多克隆抗体

中文名称：兔抗 PGBD2 多克隆抗体

英文名称： Anti-PGBD2 rabbit polyclonal antibody

储 存： 冷冻（-20℃）

抗 原： PGBD2

宿 主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The piggyBac family of proteins, found in diverse animals, are transposases related to the transposase of the canonical piggyBac transposon from the moth, <i>Trichoplusia ni</i> . This family also includes genes in several genomes, including human, that appear to have been derived from the piggyBac transposons. This gene belongs to the subfamily of piggyBac transposable element derived (PGBD) genes. The PGBD proteins appear to be novel, with no obvious relationship to other transposases, or other known protein families. The exact function of this gene is not known. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB

<b>Name of antibody:</b>	PGBD2
<b>Immunogen:</b>	Fusion protein of human PGBD2
<b>Full name:</b>	piggyBac transposable element derived 2
<b>SwissProt:</b>	Q6P3X8
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	68 kDa
<b>WB Positive control:</b>	A172 and Jurkat cell lysates
<b>WB Recommended dilution:</b>	200-1000

