

兔抗 PGBD5 多克隆抗体

中文名称：兔抗 PGBD5 多克隆抗体

英文名称： Anti-PGBD5 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： PGBD5

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The piggyBac family of proteins, found in diverse animals, are transposases related to the transposase of the canonical piggyBac transposon from the moth, <i>Trichoplusia ni</i> . This family also includes genes in several genomes, including human, that appear to have been derived from the piggyBac transposons. This gene belongs to the subfamily of piggyBac transposable element derived (PGBD) genes. The PGBD proteins appear to be novel, with no obvious relationship to other transposases, or other known protein families.
Applications:	ELISA, IHC
Name of antibody:	PGBD5
Immunogen:	Synthetic peptide of human PGBD5

Full name:	piggyBac transposable element derived 5
SwissProt:	Q8N414
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain and Human thyroid cancer
IHC Recommend dilution:	25-100

