

## 兔抗 WFDC3 多克隆抗体

中文名称：兔抗 WFDC3 多克隆抗体

英文名称： Anti-WFDC3 rabbit polyclonal antibody

别名： WAP four-disulfide core domain 3; WAP14; dJ447F3.3

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： WFDC3

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the WAP-type four-disulfide core (WFDC) domain family. The WFDC domain, or WAP signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor. The encoded protein contains four WFDC domains. Most WFDC genes are localized to chromosome 20q12-q13 in two clusters: centromeric and telomeric. This gene belongs to the telomeric cluster. Alternatively spliced transcript variants have been observed but their full-length nature has not been determined.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	WFDC3

<b>Immunogen:</b>	Synthetic peptide of human WFDC3
<b>Full name:</b>	WAP four-disulfide core domain 3
<b>Synonyms:</b>	WAP14; dJ447F3.3
<b>SwissProt:</b>	Q8IUB2
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human breast cancer
<b>IHC Recommend dilution:</b>	50-100



